

Our product and service offering

The Sonova Group is a leading global provider of hearing solutions and committed to offering the most comprehensive range of industry-leading solutions. Our product brands – Phonak, Unitron, Hansaton and Advanced Bionics – create compelling new products to offer the optimal solution for every user, provided by hearing care professionals through a broad range of service channels, including our well established global audiological care business.

Hearing Instruments Business


PHONAK

unitron

HANSATON

Cochlear Implants Business

AB Advanced Bionics



Audiological Care Business – Professional audiological services



Phonak

Audéo™ Marvel

The Phonak Audéo Marvel Receiver-In-Canal hearing aid delivers top-rated streamed sound quality¹ and automatically adapts to different listening situations thanks to AutoSense OS™. This multifunctional hearing aid is the world's first to feature technology capable of supporting binaural direct streaming from iOS, Android and Bluetooth enabled phones as well as from billions of other devices, all in stereo sound. Hearing aid wearers can enjoy the benefits of rechargeable technology and a host of dedicated apps which allow them to experience remote hearing adjustments and real-time voice-to-text call transcription.

- 1) Legarth, S., Latzel, M. & Rodrigues, T. (2018). Media streaming: The sound quality wearers prefer, Phonak Field Study News, retrieved from www.phonakpro.com/evidence, accessed October 17, 2018.

Android, Google Play and the Google-Play-Logo are trademarks of Google Inc.

IOS is a trademark of Cisco, Inc.

The Bluetooth® word mark and logos are registered trademarks owned by the Bluetooth SIG, Inc.



Unitron

Moxi™ Jump R

Moxi Jump R – Unitron’s first lithium-ion rechargeable Receiver-In-Canal hearing aid has impeccable natural sound and speech clarity. It offers a full day of hearing, including phone and media streaming, to both ears – hands-free. It’s easy to use, with no battery door – and, it looks great. With the revolutionary Unitron Flex™ ecosystem of technologies, services, and programs it is designed to make the process of buying and using a hearing aid easy, personal, and empowering.



Hansaton

sound SHD stream

The HANSATON sound SHD stream is a Receiver-In-Canal hearing system (RIC) with binaural features for a more natural hearing experience. Safe hands-free phone calls, for example in the car or while cooking, are possible for all operating systems, in short MFA (“Made For All”), for various manufacturers of Bluetooth® capable telephones. The wearers are offered an all-rounder hearing system that ideally meets the needs and demands of the digital communication age.

- 1) Smartphones and traditional mobile phones with Bluetooth® 4.2 wireless technology and most older Bluetooth® versions, Bluetooth® is a trademark owned by the Bluetooth SIG



Advanced Bionics

HiRes™ Ultra 3D Cochlear Implant

Advanced Bionics this year introduced the new breakthrough HiRes™ Ultra 3D implant, featuring a new magnet design. More and more people undergo an MRI scan (magnetic resonance imaging) during their lifetime, but this has long been a hassle as the powerful magnetic field generated by the scanner caused a strong force on the magnet contained in traditional implants. This resulted in the need to surgically remove the magnet in order to ensure a pain-free scan and meant that recipients were not able to hear until it was replaced and the healing process was complete. With the HiRes Ultra 3D implant, recipients can now undergo an MRI scan without feeling pain or discomfort and without the need to remove the magnet, thereby improving the quality of life for recipients.

